




mylife CamAPS FX app on your iPhone

iOS specific information and instructions

The Automated Insulin Delivery (AID) system mylife Loop consisting of the mylife CamAPS FX app, the mylife YpsoPump insulin pump and the FreeStyle Libre 3 / FreeStyle Libre 3 Plus or Dexcom G6 CGM (Continuous Glucose Monitoring) sensor is **now available on iPhones!**

The smartphone-based therapy management is the future of diabetes treatment. That is why we are extending the freedom of choice, as people with type 1 diabetes can choose between the iOS or Android operating systems to use the mylife CamAPS FX app. Previously available only for Android phones, the app is now accessible for iPhone users, opening the door for even more people to benefit from this technology – entirely via their smartphone.



How to download the iOS version of the mylife CamAPS FX app

Download the mylife CamAPS FX app (mmol/L or mg/dL) from the Apple App Store to your iPhone. Use the QR code to directly download the app. It is recommended to disable automatic updates in the settings for your App Store to ensure that the mylife CamAPS FX app works smoothly at all times.



mmol/L



mg/dL



Compatibility note

Before downloading the mylife CamAPS FX app, check if your smartphone with its operating system is compatible with the app.

- For minimum operating system requirements of the mylife CamAPS FX app please see <https://camdiab.com/faq>
- Smartphones compatible with FreeStyle Libre 3 / FreeStyle Libre 3 Plus: www.diabetescare.abbott/support/manuals.html ("Mobile Device & OS Compatibility" guide).
- If using a FreeStyle Libre 3 or FreeStyle Libre 3 Plus sensor, please ensure you have a smartphone with NFC (near field communication).
- Smartphones compatible with Dexcom G6 CGM system: <https://www.dexcom.com/compatibility/g6>.



More freedom. More confidence. With mylife.

How to switch between Android and iOS

If you want to switch from the Android version of the mylife CamAPS FX app to the iOS version (or vice versa), you have to keep the following information in mind.



Plan your switch carefully

Plan your switch carefully to ensure all important information is secured, the transition is as smooth as possible, and your therapy will be not interrupted. It's recommended to make the switch together with your diabetes care team. For any questions or support, feel free to reach out to your local customer service.



Backup your data

To avoid any data loss, generate and save a PDF report of your historical data before switching between the Android and iOS app. This will help you and your diabetes team maintain a clear overview of your settings like bolus calculator settings, total daily dose and personal glucose target. Use the generated PDF report to complete the setup of the mylife CamAPS FX app on your new smartphone.



Re-learning of the algorithm

You can log in with the same mylife CamAPS FX account for iPhones and Android smartphones. The mylife CamAPS FX algorithm will need to re-learn your individual insulin requirements and patterns in your data when switching between Android and iOS. This process may take a few days to adapt the insulin delivery fully to your needs.



Delete the mylife YpsoPump and the Dexcom G6 transmitter in the Bluetooth® settings of your Android smartphone, before setting up the mylife CamAPS FX app on your iPhone (or vice versa). Delete the mylife CamAPS FX app from your old smartphone.


Alert settings

- Time-Sensitive alerts: must be turned ON for the mylife CamAPS FX app to ensure alerts are received
 - Settings → CamAPS FX → Notifications → Time-Sensitive notifications
- Silent mode: using the side button will **mute ALL** alerts **including Urgent Low**.
- Focus mode/Do not disturb: will **mute ALL** alerts **including Urgent Low** unless the mylife CamAPS FX app is excluded.
- If no personal alert settings are entered, the app will automatically apply default values for all alerts and the personal glucose target.


App updates

- Automatic iOS and app updates: this may prevent the app from working until the iPhone is unlocked and/or the app reopened.
- Turn off automatic iOS updates
 - Settings → Software Update → Automatic Updates
- Turn off automatic app updates
 - Settings → Apps → App Store → App Updates
 - Note: this disables auto-updates for all apps, not just mylife CamAPS FX.
- User control: users will receive advance warning before updates are released so you can choose when to update the app.

Screen overview



Bolus calculator of the mylife CamAPS FX app on iPhones



1. Tap on the knife and fork symbol at the bottom middle of the home screen.
2. Enter the carbohydrates. Check the suggested bolus and tap “Next”.
3. Tap “Deliver” to send the bolus to the pump.
4. A 5 second count-down occurs before the bolus delivers. You can cancel it by tapping “Cancel”.

Disclaimer: The product images are for illustrative purposes only. // mylife and Ypsomed are registered trademarks of Ypsomed AG in several countries. // CamAPS is a registered trademark of CamDiab Ltd. // Dexcom and Dexcom G6 are registered trademarks of Dexcom, Inc. in the United States and/or other countries. // The sensor housing, FreeStyle, Libre, and related brand marks are marks of Abbott and used with permission. // Apple and iPhone are trademarks of Apple Inc., registered in the US and other countries and regions. // The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc., and any use of such marks by Ypsomed is under license. // Other trademarks are the property of their respective owners. // mylife Loop on iOS: The innovation described is available in selected countries and expanding to further countries soon. Expansion is contingent upon local regulatory approval.

For safety information on mentioned products, see mylife-diabetescare.com/safety