



Ypsomed Infusion Set

Fill Rate Guide

YpsoPump® Orbit®micro Infusion Sets

High-quality 31G stainless steel cannula for great comfort and reliability.



YpsoPump® Orbit®micro Infusion Sets

Available cannula lengths	Fill volume
YpsoPump® Orbit®micro 5.5 mm cannula	0.2 units (0.002 ml)
YpsoPump® Orbit®micro 8.5 mm cannula	0.2 units (0.002 ml)

Available tubing lengths	Fill volume
YpsoPump® Orbit®micro 45 cm (18") tubing	10 units (0.10 ml)
YpsoPump® Orbit®micro 60 cm (24") tubing	12 units (0.12 ml)
YpsoPump® Orbit®micro 80 cm (31") tubing	15 units (0.15 ml)
YpsoPump® Orbit®micro 110 cm (43") tubing	18 units (0.18 ml)

YpsoPump® Orbit®soft Infusion Sets

Flexible soft cannula with a low penetration force for insertion comfort.



YpsoPump® Orbit®soft Infusion Sets

Available cannula lengths	Fill volume
YpsoPump® Orbit®soft 6 mm cannula	0.3 units (0.003 ml)
YpsoPump® Orbit®soft 9 mm cannula	0.4 units (0.004 ml)

Available tubing lengths	Fill volume
YpsoPump® Orbit®soft 45 cm (18") tubing	10 units (0.10 ml)
YpsoPump® Orbit®soft 60 cm (24") tubing	12 units (0.12 ml)
YpsoPump® Orbit®soft 80 cm (31") tubing	15 units (0.15 ml)
YpsoPump® Orbit®soft 110 cm (43") tubing	18 units (0.18 ml)

YpsoPump® Inset Infusion Sets

All-in-one concept – the Infusion set with the integrated inserter and a flexible soft cannula.



YpsoPump® Inset Infusion Sets soft cannula

Available cannula lengths	Fill volume
YpsoPump® Inset 6 mm soft cannula	0.1 units (0.001 ml)
YpsoPump® Inset 9 mm soft cannula	0.2 units (0.002 ml)

Available tubing lengths	Fill volume
YpsoPump® Inset 46 cm (18") tubing	10 units (0.10 ml)
YpsoPump® Inset 60 cm (24") tubing	12 units (0.12 ml)
YpsoPump® Inset 80 cm (31") tubing	14 units (0.14 ml)
YpsoPump® Inset 110 cm (43") tubing	18 units (0.18 ml)



Orbit®

YpsoPump® Orbit®micro

Infusion set with a very thin (31G) stainless steel cannula for a 90° insertion

- Option for individuals with a history of occlusions due to bent soft cannulas and those who have an allergy to soft cannulas
- Freely rotates 360° – for more convenience and for flexibility in the infusion site selection
- With the YpsoPump® adapter – for an independent cartridge and infusion set change
- High-quality steel cannula for great comfort and reliability.
- Tan and semi-transparent tape – for more discreteness.
- Infusion set change can be performed easily and quickly
- Reusable insertion device available – the Orbit®Inserter

YpsoPump® Orbit®soft

Infusion set with a flexible soft cannula for a 90° insertion

- Freely rotates 360° – for more convenience and for flexibility in the infusion site selection
- With the YpsoPump® adapter – for an independent cartridge and infusion set change
- Tan and semi-transparent tape – for more discreteness
- Skin-friendly and breathable tape – for more comfort
- Infusion set change can be performed easily and quickly
- Reusable insertion device available – the Orbit®Inserter

YpsoPump® Inset

Infusion set with a flexible soft cannula for a 90° insertion

- All-in-one concept – the infusion set with an integrated insertion device
- Active needle safety system – for your safety
- With the YpsoPump adapter – for an independent cartridge and infusion set change
- Skin-friendly adhesive – for more comfort
- Disconnection at site – for increased flexibility
- Visible inspection window to see insertion point
- Available in grey